

Dopamine ELISA kit – Ultra-Sensitive

Any sample

Ref: BA-E-5300R

This research Dopamine ELISA kit enables the ultra-sensitive measurement of Dopamine (DA) in any biological sample. Highly flexible, the assay features a minimal sample volume as low as 1 µl, with a limit of detection at 3.3 pg/ml.

Sample type	1 - 750µl, all biological fluids
Samples/kit	40 samples in duplicate
Sensitivity	3.3pg/ml
Range	0.5 - 80ng/ml
Assay time	Extraction, acylation and enzymatic conversion (4 hours), ELISA overnight
Reactivity	Reacts with all species
Citations	Cited in 30+ papers

INFORMATIONS

Product overview

Product name	Dopamine ELISA kit – Research
Description	Enzyme immunoassay (ELISA) for the ultra-sensitive determination of Dopamine in biological samples. Research use only
Format	96-well plate
Samples	Any sample
Sample volume	1 – 750µl
Reactivity	Reacts with all species
Standard range	0 / 0.5 – 80ng/ml
Sensitivity	3.3pg/ml
Specificity	No significant cross-reactivity was observed with Dopamine analogs.
Assay time	Sample preparation 4 hours and ELISA overnight
Storage	Store at 2-8°C for to 6 months
Datasheets	Instructions for use , Material safety datasheet

PROTOCOLS

Sample preparation

Sample extraction and acylation (1h45), enzymatic conversion (2h)

ELISA

Dopamine antiserum overnight incubation, revelation and read steps (1h).

Detailed protocol

[Download instructions for use](#)

REFERENCES

Product citations

Product pictures



Dopamine ELISA kit



Dopamine ELISA kit



Dopamine ELISA kit

Contact information

Immusmol
229 Cours de l'Argonne
33 000 Bordeaux - France
Tel: +33 (0) 5 6431 1170
www.immusmol.com

To order, review, ask for technical support, visit product page at:

<https://www.immusmol.com/shop/dopamine-elisa-kit-any-sample/>