

AA-E-1100

Fast DHEA-S ELISA kit - Human serum

This Fast DHEA-S ELISA kit enables the quantitative measurement of steroid hormone Dehydroepiandrosterone sulfate in human serum samples within 60 minutes, to evaluate reproductive health and adrenal function. Intended for *in vitro* diagnostics and research applications, the test comes in 96-well plate format, allowing the screening of large experimental series within short lead times.



Sample type	Serum
Capacity	96-well plate
Sensitivity	0.005 µg/mL
Range	0/0.005 – 10 µg/mL
Assay time	60 min
Reactivity	Human
Application	Research use or in vitro diagnostics (IVD in EU only)
Storage	Store at 2-8°C for to 6 months
Sample collection & storage	Serum: Store samples at 2-8°C for up to 48h or -20°C for longer period (up to 6 months).
Sample preparation	Enzymatic conjugation (45 min)
ELISA	Incubation, revelation and read steps (15-20 min)

To see data, publications, and download instructions for use and MSDS, visit:
<https://www.immusmol.com/product/fast-dhea-s-elisa-kit-serum/>