

FC-E-3100

## Food Safety Histamine ELISA kit - Any food or beverage

This fast Food Safety Histamine ELISA kit is tailored for the measurement of histamine in food and beverage samples within 75 minutes (including fish and fish products, sausages, dairy products and wine). Simple and rapid, the test displays high correlation with HPLC. **This product is only available for customers outside of France.**



<b>Sample type</b>	1g fish meal; 10g fresh fish, sausage or cheese; 10 µl milk; 20 µl wine
<b>Capacity</b>	96-well plate
<b>Sensitivity</b>	0.15µg/l
<b>Range</b>	Depends on type of sample
<b>Assay time</b>	75 min
<b>Reactivity</b>	Reacts with all species
<b>Application</b>	Food safety
<b>Storage</b>	Store at 2-8°C for to 6 months
<b>Sample collection &amp; storage</b>	See instructions for use
<b>Sample preparation</b>	Sample acylation (15 min)
<b>ELISA</b>	ELISA incubation, revelation and read steps (1h)

To see data, publications, and download instructions for use and MSDS, visit:  
<https://www.immusmol.com/product/food-safety-histamine-elisa-kit/>