

BA-E-2500

GABA ELISA kit - Any sample

This GABA ELISA kit enables the measurement of Gamma-aminobutyric acid in any biological sample. Intended for *in vitro* diagnostics, this quantitative assay can be performed in cell cultures, tissue homogenates, urine, plasma, serum, etc.

[Cited in 18 papers](#)



Sample type	100 – 300µL, all biological samples
Capacity	96-well plate
Sensitivity	Depending on sample matrix and volume
Range	25 – 2 500 ng/mL 75 – 7 500 ng/mL
Assay time	Extraction & derivatization (2h30), incubation overnight, ELISA revelation & read step (1h)
Reactivity	Reacts with all species
Application	In vitro diagnostics (in EU-only), or Research use
Storage	Store at 2-8°C for to 6 months
Sample collection & storage	See instructions for use
Sample preparation	Sample extraction and derivatization (2h30)
ELISA	GABA antiserum overnight incubation, revelation and read steps (1h)

To see data, publications, and download instructions for use and MSDS, visit:
<https://www.immusmol.com/product/gaba-elisa-kit-any-sample/>