

BA-E-1100R

Histamine release in whole blood - ELISA kit Add-on (XL)

Intended for research use, this add-on to the [BA-E-1000R histamine ELISA kit](#) enables the measurement of histamine release in heparinized whole blood samples, allowing for instance to determine the sensitivity of basophilic granulocytes to specific allergens. Important: this product is only sold in 10-unit format.



Sample type	Heparinized whole blood
Capacity	10 x 96 determinations
Range	0/0.5 – 50 ng/mL
Assay time	Sample preparation 90 min
Reactivity	Reacts with all species
Application	Research use only
Storage	Store at 2-8°C for to 6 months
Sample collection & storage	See instructions for use

To see data, publications, and download instructions for use and MSDS, visit:
<https://www.immusmol.com/product/histamine-release-kit-whole-blood-samples/>