

IS1030

L-DOPA Antibody - Rabbit Polyclonal

The anti-L-DOPA antibody IS1030 was validated for immunofluorescence in crayfish tissues and mouse primary neurons using the [STAINperfect immunostaining kit A](#), demonstrating specific and accurate labeling of L-DOPA.



Clonality	Polyclonal antibody
Host	Rabbit
Reactivity	Reacts with all species
Tested samples	Mouse primary neurons and crayfish tissues
Staining procedure	STAINperfect immunostaining kit A (protocols for cell cultures, whole mounts, tissue sections)
Format	50µL (approx. 40 tissue sections)
Storage	Store at +4°C for short term (1-2 months). Aliquot and store at -20°C for long term. Avoid repeated freeze / thaw cycles
Immunofluorescence (IF)	Dilute antibody with the antibody diluent provided in the STAINperfect immunostaining kit A . Use at 1/250 -1/1000 dilution. Follow the STAINperfect protocol suited to your sample (Cell cultures, Whole mounts, Tissue sections)
Comments	Optimal working dilutions must be determined by the end-user
Restrictions	For research use only

To see data, publications, and download instructions for use and MSDS, visit:
<https://www.immusmol.com/product/l-dopa-antibody-rabbit-pab/>