

IS1076

L-Leucine Antibody - Rabbit Polyclonal

The anti-L-Leucine antibody IS1076 enables the detection of the essential amino acid L-Leucine in biological samples. Our team has validated the affinity and specificity of this polyclonal antibody. Compatible with [STAINperfect immunostaining kit A](#), the antibody can be used for IF/ICC/IHC imaging applications.



Clonality	Polyclonal antibody
Host	Rabbit
Reactivity	Reacts with all species
Eligible samples	Cell cultures, whole mounts, tissue sections
Staining procedure	STAINperfect immunostaining kit A
Format	50µL (approx. 40 tissue sections)
Storage	Store at +4°C for short term (1-2 months). Aliquot and store at -20°C for long term. Avoid repeated freeze / thaw cycles
Immunofluorescence (IF)	Dilute antibody with the antibody diluent provided in the STAINperfect immunostaining kit A . Use at 1/250 -1/1000 dilution. Follow the STAINperfect protocol suited to your sample (Cell cultures, Whole mounts, Tissue sections)
Comments	Optimal working dilutions must be determined by the end-user
Restrictions	For research use only

To see data, publications, and download instructions for use and MSDS, visit:
<https://www.immusmol.com/product/l-leucine-antibody-rabbit-polyclonal/>