

BA-E-5800R

Ultra-Sensitive Histamine ELISA kit - Any sample

This Ultra-Sensitive Histamine ELISA kit enables the accurate measurement of histamine in any biological sample: tissue homogenates, cell culture supernatants, body fluids, etc. With a minimal sample volume as low as 10µL, this histamine immunoassay allows longitudinal studies in mice and rat models.

[Cited in 7 papers](#)



Sample type	Any sample
Capacity	96-well plate
Sensitivity	0.2 ng/mL
Range	0.5 – 50 ng/mL
Assay time	Sample acylation (90 min), incubation (3h or overnight), ELISA revelation & read steps (90 min)
Reactivity	Reacts with all species
Application	Research use only
Storage	Store at 2-8°C for to 6 months
Sample collection & storage	See instructions for use
Sample preparation	Sample predilution and acylation (1h30)
ELISA	Histamine antiserum incubation (3h or overnight), revelation and read steps (1h30)

To see data, publications, and download instructions for use and MSDS, visit:
<https://www.immusmol.com/product/ultra-sensitive-histamine-elisa-kit/>