

BA-E-5900R

Ultra-Sensitive Serotonin ELISA kit, Any sample

This Ultra-Sensitive Serotonin ELISA kit enables the accurate measurement of 5-hydroxytryptamine (5-HT) in any biological sample. Intended for research use only, the assay features a lower limit of detection at 5pg/ml and a minimal sample volume as low as 1µl.

[Cited in 27 papers](#)



Sample type	Any sample
Capacity	96-well plate
Sensitivity	5pg/mL
Range	0.015 - 2.5ng/mL
Assay time	Sample acetylation (30 min), overnight incubation, ELISA revelation & read step (1h)
Reactivity	Reacts with all species
Application	Research use only
Storage	Store at 2-8°C for to 6 months
Sample preparation	Sample acylation (30min)
ELISA	Antisera overnight incubation, revelation and read steps (1h)

To see data, publications, and download instructions for use and MSDS, visit:
<https://www.immusmol.com/product/ultra-sensitive-serotonin-elisa-kit/>